

# Careturmer<sup>®</sup>

## INNOVATIVE CARE TECHNOLOGY

### 3 in 1 solution

- Carefully turning patients
- Preventing development of pressure ulcers
- Optimizing working environment

Careturmer by GDV Technology




## Careturmer is turning people with care

Careturmer is a unique and innovative medical device designed to reposition patients automatically and to support care staff with daily mobilization.

Careturmer turns the patient gently from side to side throughout the night and provides pressure distribution to prevent development of pressure ulcers.

Careturmer is optimized for single caregiver working environments. It optimizes the working environment and saves resources.



### Careturmer fits following medical beds:

- Stiegmeyer Elvido-Movo
- wissner-bosserhoff Sentida
- Invacare SB 755
- Lojer ScanAfia XHS

More versions are in process and can be made on demand.



## Careturner turns the patient gently from side to side

Careturner is easy to install on a medical bed, and works in conjunction with the beds existing movements.

- **Manual mode:** With a push on a button it turns the enduser from one side to the other to support the care staff in daily moving and handling procedures. In manual mode the wings can elevate up to 80°.

- **Automatic mode:** The control box has a preprogrammed setting and option for two additional programs. Careturner operates very silent and without any vibration allowing the patient to sleep undisturbed. In automatic mode the wings can elevate up to 30°.

### Intended use of Careturner

- Patients at nursing homes, hospitals, hospices or at home
- Patients suffering from dementia
- Bariatric patients
- Patients in intensive care units
- Chronic disabled patients suffering from e.g. rheumatism, sclerosis or brain injury

## Careturner is tested and approved by TÜV with patient weight up to 165 kg



Careturner works in conjunction with the beds existing movements



The programmable control box



Sides can be angled up to 60-80°. Sling can be used by 1 caregiver



"Hug & Cradle" function gives peace and calms down the patient

## About GDV Technology A/S

GDV Technology A/S is a Danish manufacturer of medical equipment for immobile patients. The company was founded in 2014.

Careturner is a patented solution developed in collaboration with Danish doctors, nurses and experts. The product was first introduced to Danish nursing homes in 2016. Careturner was tested and approved by TÜV in 2017.

Careturner was elected start-up of the year in the category Care at Altenpflege, Germany in 2018 and has also received prizes at the Occupational Therapy Show, UK 2018 and Health & Rehab, Scandinavia in 2016.



Scan QR-code and watch product demonstration video



GDV Technology A/S  
Værkstedsgården 15, st. th.  
DK-2620 Albertslund  
Tel. +45 2097 8115  
Mail: info@gdvtech.com  
www.careturner.com